



UV-B Narrowband TL

TL 40W/01 RS

More than 400 independent clinical studies have proven that the UVB Narrowband TL lamps are safer and more effective than any other lamps in their class. That is because these lamps emit only a very narrow waveband from the 'B' bandwidth of the UV spectrum (290 to 315). This narrow waveband is between 305 and 315 nm and peaks at 311 nm: the most efficacious waveband for the treatment of psoriasis. This means that treatment is much more focused and exposure times are much shorter. This in turn leads to a reduction of side effects such as reddening of the skin and itching. All of this makes them ideal for phototherapy treatment of diseases such as psoriasis and vitiligo. What's more, because the overall dosage of this narrowband radiation can be closely controlled, these lamps are suitable for home therapy. N.B.: Our UVB lamps are NOT registered with FDA as medical devices as they are NOT packaged or labeled for commercial distribution for health-related purposes. US customers are referred to the UVB and UVA lamp range brochure US version.

Product data

• General Information

Cap-Base	G13 [Medium Bi-Pin Fluorescent]
Bulb Shape	T38 [T 38mm]
Main Application	Medical Therapy
Life To 50% Failures (Nom)	9000 h
Useful Life (Nom)	9000 h
System Description	Rapid Start

• Light Technical

Color Code	01
Color Designation	Ultra Violet B
Chromaticity Coordinate X (Nom)	208
Chromaticity Coordinate Y (Nom)	192
Lumen Depreciation At 500 Hours	10 %
Lumen Depreciation At 1000 Hours	15 %

• Operating and Electrical

Power (Rated) (Nom)	39 W
Lamp Current (Nom)	0.43 A
Voltage (Nom)	104 V

• UV

UV-B Radiation 100 hr (IEC)	7.5 W
UV-B Radiation 5hr (IEC)	7.7 W

• Product Data

Full product code	871869666239700
Order product name	TL 40W/01 RS
Order code	928011300130
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	25
Material Nr. (12NC)	928011300130

PHILIPS

UV-B Narrowband TL

Net Weight (Piece) 292.000 g

Warnings and Safety


- A lamp breaking is extremely unlikely to have any impact on your health. If a lamp breaks, ventilate the room for 30 minutes and

remove the parts, preferably with gloves. Put them in a sealed plastic bag and take it to your local waste facilities for recycling. Do not use a vacuum cleaner.

Dimensional drawing



TL 40W/01 RS

Product	D	A	B	B	C
TL 40W/01 RS	40.5 mm	1199.4 mm	1206.5 mm	1204.1 mm	1213.6 mm



TL G13

Photometric data



© 2016 Philips Lighting Holding B.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. (Royal Philips) or their respective owners.

www.philips.com/lighting

2016, October 11
data subject to change